



N'DURA HOOF™ INSTRUCTIONS FOR USE
PROSHIELD™ AND COLORSHIELD™

PRECAUTIONS

- Precisely follow directions.
- Store in a cool place.
- Avoid all skin contact. If N'Dura Hoof contacts skin, wash off with soap and water.
- Keep tightly sealed.
- Do not expose applicator or can contents to light or hardening will occur.
- Keep out of reach of children.
- Wipe can threads with alcohol after every use to prevent sticking.
- Remove uncured N'Dura Hoof with acetone.

INSTRUCTIONS FOR USE

1. Apply N'Dura Hoof products in the shade only.
2. Apply to one hoof at a time.
3. Smooth, sand or rasp hooves and brush off any loose particles.
4. Make sure hoof is clean, dry and free of any oil or grease.
5. If hoof is oily, wipe with isopropyl alcohol. **DO NOT USE ACETONE.**
6. Dry the hoof with external heat if necessary.
7. Apply N'Dura Hoof in a thin layer to the exterior hoof wall, avoiding the coronet band.
8. Smooth application as desired.
9. Expose N'Dura Hoof coating to the N'Dura Hoof UV LED Light:
 - a. ProShield for 30 seconds
 - b. ColorShield for 60 seconds
10. OPTIONAL: Repeat, if thicker coating or more color is desired.

REMOVAL

1. Sand or file off at any time.

INGREDIENTS: Proprietary monomers, photo initiators, thickeners, pigments and glitter.

US PATENTS PENDING

www.NDuraHoof.com

1.858.348.1122